

A.M.I.® OptiClean Device

Product Group Laparoscopy | Issue 02/2008



> **Minimise risk by maximising vision**

- Cleans the lens at the touch of a finger
- Saves up to 20 minutes operating time
- Makes scope removal unnecessary
- Gives you a clear view and complete control

A.M.I.® OptiClean Device

A.M.I. OptiClean - a clear view in seconds!

- Prevents the lens fogging up during endoscopic surgery
- Enables fast lens cleaning without scope removal
- Attaches easily to 10 mm laparoscopic scopes
- Saves up to 20 minutes during advanced procedures
- Allows unobstructed vision at all times

| Order Code | A.M.I. OptiClean Device | Technical Details |
|------------|--|--|
| AOC 400 | A.M.I. OptiClean Device 305 consisting of: <ul style="list-style-type: none">- A.M.I. OptiClean snap-on profile- Tubing- Irrigation syringe- Suction tube puncture needle- Suction tube lock | Suitable for 10 mm 0° to 30° scopes A 10 mm scope with OptiClean attached requires a 12 mm trocar Profile length: 305 mm Single-use, sterile-packed |
| AOC 402 | A.M.I. OptiClean Device 285 consisting of: <ul style="list-style-type: none">- A.M.I. OptiClean snap-on profile- Tubing- Irrigation syringe- Suction tube puncture needle- Suction tube lock | Suitable for 10 mm 0° to 30° scopes A 10 mm scope with OptiClean attached requires a 12 mm trocar Profile length: 285 mm Single-use, sterile-packed |

A.M.I. Headquarters:

A.M.I. GmbH
Im Letten 1
6800 Feldkirch . Austria
t +43 5522 90505-0
f +43 5522 90505-4006
e info@ami.at
www.ami.at